



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET MEDAP P7030 surgical suction (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Manufactured in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 88 W,
MAQUET MEDAP 7030 surgical suction device for aspirating secretions and fluids during procedures. Primarily used in outpatient clinics and hospital emergency departments (ED),
Specification:
Generated vacuum: -76 kPa/ -675 mmHg,
Maximum throughput: 12 l/min,
Smooth adjustable suction strength,
Built-in manometer,
Suction device mounted on a mobile dedicated trolley for easy transport,
Included:
Autoclavable jar (121°C) SERRES with a capacity of 2 L and disposable insert,
Suction hose length: 150 cm,
Power cord,
Overall dimensions: 40 x 30 x 85 cm,
Weight: 9 KG,
Has current inspection and is ready for use,
Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of warranty extension up to 12 months for both domestic and international markets, with surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

