



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LIEBHERR FKS 1802 refrigerated cabinet (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Manufactured in Austria,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 130 W,
Specification:
Silent cooling,
Mechanical control,
Temperature from +1°C to +7°C,
Exterior: white steel sheet,
Interior: white plastic,
2 indents for transport handle,
Glass doors that open,
Included:
3x shelves with dimensions: 41 x 54 cm,
Dimensions:
External: 60 x 60 x 88.5 cm,
Internal: 51.3 x 44.1 x 73.4 cm,
Weight: 45 KG,
Has a current inspection and is ready for operation,
Technical Passport (Service Report) issued and valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Option to extend warranty for an additional charge up to 12 months for domestic and international markets,
Financing option (Poland only): Installments, Leasing, Loan,

