



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ZEISS IOL MASTER optical biometer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 90 VA,
- Software: 5.4.,
- **Characteristics:**
 - Fast and precise measurements of axial length, corneal radii, anterior chamber depth, and white-to-white in a single session,
 - All measurements are non-invasive, ensuring maximum patient comfort,
 - It uses the patented partial coherence interferometry (PCI) method to measure axial length,
 - Intraocular lens power calculations based on internationally recognized formulas: Haigis, HofferQ, Holladay, SRK II and SRK/T,
 - Possibility of optimizing fixed lenses based on postoperative refraction results,
 - IOL database with the ability to add data on selected lenses,
- **Technical specifications:**
 - Axial length:
 - Area: 14 – 40 mm,
 - Display scaling: 0.01 mm,
 - Keratometer:
 - Area: 5 – 10 mm,
 - Display scaling: 0.01mm,
 - Anterior chamber depth:
 - Area: 1.5 – 6.5 mm,
 - Display scaling: 0.01 mm,
 - White from white:
 - Area: 8 – 16 mm,
 - Display scaling: 0.1 mm,
 - Chin rest adjustment range: 20 – 26 cm
- **Possibility to select menu in: GER, ENG, ESP, FR, IT, SWE,**
- **The equipment is placed on a mobile electric table with the possibility of adjusting the tabletop height in the range of: 70 - 100 cm,**
- **Included in the set:**
 - CARL ZEISS mobile electric table,
 - IOL MASTER artificial eye,
 - HP LASERJET 1020 printer,
 - KEYSONIC Wired Multifunction Keyboard,
 - User manual in English (PDF),
 - Power cord,
- Dimensions:
 - Tabletop: 68 x 38 cm,
 - Biometer: 32 x 40 x 58 cm,
 - Total: 70 x 60 x 132/162 cm,
- Weight: 50 KG,
- Cat 02,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
 - **The warranty period does not cover light bulbs and other light sources installed in the device,**
- Financing options (Poland only): Installments, Leasing, Loan,

