



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical records cabinet (Reconditioned)

Refurbished (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

Medical records cabinet

is intended for safe and organized storage of patient files and other medical documents. It is equipped with sliding drawers adapted for segregation of files with built-in guides allowing smooth opening,

Solid metal construction ensures durability and protection of documents,

Features:

Freely sliding drawers,

Metal solid construction,

Built-in central lock,

In the set:

Keys for locking drawers (2 pcs),

Dimensions:

Drawers: 26 x 58 x 22 cm,

Total: 35 x 64 x 128 cm,

Weight: 55 KG,

Has an up-to-date inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend the warranty for an additional charge up to 12 months for the domestic and international market,

Financing options (Poland only): Installments, Leasing, Loan,

