



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDAP P2050 surgical suction (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 180 W,
MEDAP P2050 surgical suction is a device for suctioning secretions and fluids during procedures. It is mainly used in treatment clinics and hospital emergency departments (ED),
Specification:
Generated vacuum: 0.96 bar,
Flow rate: 35 l/min,
Smooth adjustment of suction force,
Built-in manometer,
Suction unit mounted on a mobile dedicated cart for easy transport,
In the kit:
Power cord,
Autoclavable jar with capacity of 1 L,
Overall dimensions: 55 x 38 x 80 cm,
Weight: 9 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months on the domestic market (Poland),
3 months on the international market,
Possibility to extend the warranty for an additional cost up to 12 months for both domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

