



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Ophthalmic unit MECCANOTTICA MAZZA (Reconditioned)**

Refurbished (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
Italian production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,

Ophthalmic unit Meccanottica Mazza

is a modern diagnostic workstation intended for comprehensive ophthalmic and optometric examinations,

The device integrates in one construction elements essential for the work of a doctor or optometrist,

The unit's design was created with ergonomics in mind, easy access to instruments and patient comfort. A stable base and worktop locking system ensure safety during examinations,

Specification:

Electric height adjustment of the patient chair in the range 48 - 63 cm,

Possibility of rotating the patient chair by 360° with locking option,

Rotating top in the range 90° (with electric rotation lock),

2 built-in light sources (12 V / 55 W),

Possibility to equip the unit with a selected slit lamp / auto-refractor from our offer for an additional fee,

Included in the set:

Power cable,

Dimensions:

Patient chair:

Seat: 46 x 43 cm,

Backrest: 48 x 44 cm,

Top: 94 x 40 cm,

Total: 120 x 100 x 205 cm,

Total weight:

Has current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend warranty for an extra charge to 12 months for both domestic and international markets,

Warranty period does not cover bulbs and other light sources installed in the device,

Financing options (only Poland): Installments, Leasing, Loan,

