



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **SCHMITZ ARCO MEDI-MATIC 114.495 gynecological chair (Reconditioned)**

Refurbished (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
New upholstery with medical certificate in cream color,  
ARCO gynecological chair  
is intended for gynecological examinations and procedures,  
Technical specification:  
Electrical tilt adjustment of the chair within 50° using a foot controller,  
Adjustable leg supports within 240°,  
Retractable waste tray,  
Included:  
Power cord,  
Wired foot controller,  
Retractable waste tray with dimensions: 29 x 32 cm,  
2 leg rests,  
Dimensions:  
Headboard diameter: Ø30 cm,  
2-segment backrest:  
Top: 56 x 30 cm,  
Bottom: 57 x 55 cm,  
Overall dimensions: 160 x 80 x 140 cm,  
Maximum load: 180 KG,  
Weight: 115 KG,  
Has a current inspection and is ready for use,  
Technical Passport (Service Report) issued and valid for 12 months,  
Warranty:  
6 months on the domestic market (Poland),  
3 months on the international market,  
Option to extend the warranty for an additional fee up to 12 months on the domestic and international market,  
Financing options (Poland only): Installments, Leasing, Loan,

