



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Gynaecology chair SCHMTIZ MEDI-MATIC ARCO 114.595 (Reconditioned)**

Refurbished (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Gynaecology chair ARCO  
with a single electric motor, intended for gynecological examinations and procedures,  
New upholstery with medical certificate, in green color,  
Chair features:  
Elegant design, modern technology and favorable quality-to-price ratio,  
Backrest upholstery is removable, underneath there is a paper roll chamber with a width of up to 40 cm,  
Trendelenburg position adjustment in the range: 7°,  
Ergonomic leg supports system with integrated patient handles,  
Manual height and leg supports positioning adjustment,  
Controlled by a wired remote control with the ability to memorize 3 chair positions, and a foot attachment,  
Low starting height, facilitating getting on and off the chair,  
Footrest for the doctor,  
Included in the set:  
Waste container with dimensions: 29.5 x 30 x 5 cm,  
Wired remote control,  
Single-function foot attachment (up/down),  
Dimensions:  
Headrest: Ø 30 x 6 cm,  
Seat with backrest: 87 x 56.5 x 3.5 cm,  
Maximum seat height: 132 cm,  
Overall: 155 x 90 x 140 cm,  
Maximum load: 180 KG,  
Weight: 115 KG,  
Has a current inspection and is ready for operation,  
Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months in the domestic market (Poland),  
3 months in the international market,  
Option to extend warranty for an additional fee up to 12 months for domestic and international markets,  
Financing available (Poland only): Installments, Leasing, Loan

