



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET 1130.02B Operating Table (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
.real product photos,
German production,
Power supply: 230 V,
Emergency power: via 4 new built-in 12 V / 12 Ah batteries,
Frequency: 50/60 Hz,
Power: 400 VA,
Medical-grade upholstery in black,
MAQUET Model 1130.02B
an ergonomic and stable operating table intended for all surgical fields and applications,
Controlled by a wired remote control in multiple planes:
Heights in the range: 73 - 110 cm,
Vertical tilt angle:
Trendelenburg +30°
Anti-Trendelenburg -30°,
Lateral tilt: 25°,
Backrests, split into two independent sections,
Leg segment height (removable, two-segment split),
Manual adjustment:
Headrest,
Leg segment - lateral abduction,
Features:
Base made of stainless steel,
Electrohydraulic drive,
Fluoroscopy translucent,
Arm compatible,
Rails for mounting additional accessories,
Mobile on four wheels including two steering wheels,
Manual brake,
Table shows visible signs of use not affecting work comfort,
Possibility to upgrade the table (for an extra charge) with additional accessories, example visible in the last photos

Included:
Power cord,
Wired control remote,
Dimensions:
Top surface: 204 x 54 x 5 cm,
Overall: 204 x 62 x 73/110 cm,
Maximum load: 225 KG,
Weight: 265 KG,
Has a current inspection and is ready for use,
Technical Passport (Service Report) issued and valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee up to 12 months for domestic and international markets,
Financing options (only Poland): Installments, Leasing, Loan,

