



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Operating table ALVO ETIUDA 4-02 (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Polish production,
Power supply: 230 V,
Emergency power: via 2 new built-in 12 V / 1.2 Ah batteries,
Frequency: 50/60 Hz,
Power: 350 VA,
Medical-grade upholstery in black,
5-segment,
Base made of stainless steel,
Electrical height adjustment of the chair position using a wired control panel:
Head orientation selection,
Table "0" function,
Upper-lower: 75-108 cm,
Vertical tilt - Trendelenburg / Anti-Trendelenburg position, +35° / -35°,
Lateral tilt: +20° / -20°,
STOP function - stops the current position,
Electrical adjustment using integrated panel in the table column:
Head orientation selection,
Table "0" function,
Upper-lower: 75-108 cm,
Vertical tilt - Trendelenburg / Anti-Trendelenburg position, +35° / -35°,
Lateral tilt: +20° / -20°,
STOP function - stops the current position,
Manual adjustment of settings:
Leg segment manually adjustable -100° / +20°, via gas springs, split abducted to the sides 90°,
Manual headrest adjustment +35°/-20° via gas springs,
Table top travel: 30 cm,
Manually adjustable backrest +80° / -20° via gas springs,
Manual adjustment of the lumbar section in the range: 11 cm (crank),
Manual brake,
Includes:
X-RAY transparency,
Rails for mounting additional accessories,
Overall dimensions:
Width: 50 cm,
Length: 210 cm,
Thickness: 6 cm,
Mobile table with built-in central brake,
Possibility to equip the table with additional accessories visible in the last photos,
Included in the set:
Power cord,
Wired handheld remote,
Crank for lumbar section adjustment,
Maximum load: 250 KG,
Weight: 240 KG,
Cat 03,
Has current inspection and is ready for operation,
Technical Passport (Service Report) issued valid for 12 months,
Warranty:
6 months for the national market (Poland),
3 months for the international market,
Possibility to extend the warranty for a surcharge up to 12 months for the national and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

