



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MAQUET 1131.02.B0 Operating Table (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230 V,
- Emergency power supply: via 4 new built-in 12 V / 10 Ah batteries,
- Frequency: 50/60 Hz,
- Power: 400 VA,
- **Upholstery with medical certificate in black (possibility of replacing the upholstery with a new one for an additional fee),**
- **MAQUET 1131.02.B0** is a comfortable, mobile operating table designed for all surgical, orthopedic and surgical specialties,
- **The MAQUET** 4-segment table has a wide range of segment expansion options to facilitate the performance of individual functions,
- Electrohydraulic drive,
- **Electrically adjustable using a wired remote control in many planes:**
  - Height range: 74 - 110.5 cm,
  - Trendelenburg and anti-Trendelenburg positions in the range of: 30°,
  - Lateral tilt range: 24°,
  - Backrest angle range: -45° to + 70°,
  - "0" position button (no table height change),
- **Manual adjustment:**
  - Tilt angle and headrest in two ranges,
  - Longitudinal movement of the bed in the range of: 25 cm,
  - Leg segment with adjustable tilt,
- **It has:**
  - X-RAY translucency,
  - Rails for mounting additional equipment,
  - Power cord hidden in the base
- **Equipment:**
  - Leg support segment with the possibility of free retraction **(the support has visible signs of wear that do not affect the comfort of use)**,
- **Included in the set:**
  - Wired foot controller,
  - Leg support segment,
- Dimensions:
  - Headrest: 51 x 26 cm,
  - Backrest: 51 x 52 cm,
  - Bed segment: 51 x 58 cm,
  - Segmented leg support: 54 x 50 cm,
  - Total length with segmented support: 218 cm,
- Maximum load capacity: 250 KG,
- Weight: 240 KG,
- Cat 2,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months to the international market,**
  - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

