



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BAKMED KILLED-KM1 7.10 Surgical lamp (Ex-Display)

New (exhibit),
Actual product photos,
Manufactured in P.R.C.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 14 VA,

Surgical lamp BAKMED for medical diagnostics and lighting of the treatment field provides professional medical lighting with high and procedure-appropriate parameters,

The light source in the surgical lamp consists of 7 LED bulbs grouped in one optical head mounted on a two-arm straight boom allowing free directing of light onto the treatment area,

Surgical lamp
KILLED-KM1 7.10

is intended for use as the main or additional treatment lighting in operating rooms and treatment rooms as well as medical offices, and also in veterinary facilities,

Technical specification:

Directed non-shadowing head with 7 LED light points with a life of 35,000 h,

Light intensity 35,000 lx at a distance of 1 meter,

Light patch 180 mm at a distance of 1 meter,

Color temperature 4700K white slightly cool,

Color rendering index Ra 93,

Energy-saving - power consumption only 14 VA,

Two-arm boom with joint with a reach of 140 cm,

Vertical axis rotation range on wall mount 180°,

In the set:

Control unit with mounting,

Set of mounting screws,

Dimensions:

Total length of the 2-segment mounting arm: 153 cm,

Weight: 6.5 KG,

Has a current inspection and is ready for operation,

Technical Passport (Service Report) issued and valid for 12 months,

Warranty:

6 months on the domestic market (Poland),

3 months on the international market,

Warranty extension available for an extra charge up to 12 months for both domestic and international markets,

Warranty does not cover bulbs and other light sources installed in the device,

Financing options (Poland only): Installments, Leasing, Loan,

