



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET Operating Table - ALM 5090 (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
French production,
Power supply: 230 V,
Emergency power: via 2 new built-in 12 V / 12 Ah batteries,
Frequency: 50/60 Hz,
Power: 180 VA,
MAQUET - ALM 5090
is a comfortable, mobile operating table intended for all surgical specialties,
New upholstery (memory visco) with medical certificate in black,
4-segment,
Electromechanical drive,
Electrically adjustable using a wired remote control in multiple planes:
Height range: 78 - 113 cm,
Trendelenburg and anti-Trendelenburg position in range: 30°,
Lateral tilt in range: 18°,
Backrest angle range: from -35° to 85°,
Manual adjustment:
Headrest tilt angle range: from -90° to 60°,
Divided leg segment, abducted to the sides, adjustable tilt angle,
Features:
X-RAY translucency,
Rails for mounting additional accessories,
Mobile with electric brake (lock via remote control),
Possibility to equip the table (for an additional fee) with additional accessories, examples visible in the last photos,
Included:
Charger + power cable,
Wired remote control,
Dimensions:
Headrest: 53 x 30 cm,
Backrest: 53 x 34 cm,
Seat: 53 x 46 cm,
Leg segment: 25 x 75 cm,
Total length: 195 cm,
Max load: 190 KG,
Weight: 195 KG,
Has a current inspection and is ready for work,
Technical Passport (Service Report) issued for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Warranty extension available for a fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

