



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



MAQUET ALPHAMAXX 1133.12B3 Operating Table (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Made in Germany,
Power supply: 230 V,
Emergency power: via 2 new built-in 12 V / 17 Ah batteries,
Frequency: 50/60 Hz,
Power: 400 VA,
MAQUET ALPHAMAXX
an operating table created in close cooperation with physiotherapists, surgeons and nursing staff. Facilitates daily workflow in the operating room, enhances patient safety and comfort. It can be easily modified according to the needs of any surgical discipline in just a few quick steps,
New upholstery (memory visco) with medical certification in black,
Adjustable in multiple planes using a wired remote control or via the column panel:
Height adjustment range: 60 - 105 cm,
Way length travel: 22 cm,
Zeroing - return to horizontal blade position,
Flex and Reflex - tilt the blade in half in both directions,
Lateral tilt to the right and left in the range of: 20°,
Vertical tilt - Trendelenburg and Anti-Trendelenburg position: 30°,
Backrest adjustment: -40° / +80°,
Adjustment of the knee segment backrest in the range: 90°
Possibility to move the table forward / backward using electric drive,
Electric brake,
Manual adjustment:
Headrest: -45 / +90°,
Leg section - lateral abduction: 45°,
Specification:
4-segment,
Electrohydraulic drive,
Universal top,
X - RAY translucency,
Base housing made of durable plastic,
Equipped with rails for mounting additional accessories,
Table adapted for C-arm,
Mobile, on double wheels
Included:
User manual in English (PDF),
Power cord,
Control remote,
Possibility to retrofit the table (at additional cost) with additional accessories, examples visible in the last photos,
Dimensions:
Headrest: 50 x 28 x 6 cm,
Back support: 54 x 25 x 6 cm,
Bed segment: 54 x 85 x 6 cm,
Leg backrest segment: 22 x 56 x 6 cm,
Maximum load: 450 KG (patient + accessories),
Weight: 340 KG,
Has current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months on domestic market (Poland),
3 months on international market,
Warranty extension available for an additional fee up to 12 months on domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

