



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **DEWERT 900-90 KH Standing Table (Reconditioned)**

Refurbished (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Medical-certified blue upholstery,  
Table mechanism provides smooth, electric backrest tilt adjustment from any position. It is possible to achieve both positive and negative tilt angles. Transition time from horizontal to vertical position is optimized for oblique assessment recommendations,  
Specification:  
Electric backrest tilt adjustment via electric foot switch:  
Trendelenburg position: 8°,  
Ability to set the backrest to vertical position: 90°,  
Backrest tilt adjustment range: 34°,  
Straps with adjustable length (2 pcs),  
Rolling base on 4 braking wheels,  
Dimensions:  
Backrest: 73 x 180 cm,  
Leg rest: 50 x 38 cm,  
Backrest height in lying position: 70 cm,  
Overall: 85 x 195 x 110 cm,  
Weight: 100 KG,  
Has an up-to-date inspection and is ready for work,  
Technical Passport (Service Report) issued and valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility to extend the warranty for an additional 12 months for domestic and international markets,  
Financing options available (Poland only): Installments, Leasing, Loan,

