



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Gurney (Ex-Display)

- Condition: New: Ex-Display,
- Made in Poland,
- Very light-frame made from aluminium,
- Dimensions: W 55 cm x L 180 cm,
- Total dimensions: L 190 cm x W 55 cm x H 15 cm,
- Handles length: 19 cm (hidden in the frame),
- Pillow pocket: 37 x 35 cm.
- Warranty: 30 days

