



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Accessories (New)

### Video Camera:



The central focus of the camera in the body of the lamp ensures a largely undistorted recording of the procedure and ensures the best possible quality shadow-free film.

The EMAICAM OR camera fits perfectly into the 560 and 500 series surgical lights.

### HD Camera

#### Technical Specifications:

Image - 1/3" type CMOS,  
Resolution - 2 000 000 pix,  
Zoom (digital) - 12 x,  
Zoom (optic) - 10 x.

#### Additional Equipment:

- Ceiling-mounted set for EMALED 300,
- Ceiling-mounted set for the models 500/600,
- Frame for indirect suspension,
- Power cord: 230 V,
- Wireless controller,
- Wall-mounted control panel,
- Uninterruptible power supply with storage battery,
- Battery for the ceiling-mounted lamps 500 i 560,
- Charger for ceiling-mounted lamps 500 i 560,
- Camera SONY SD with remote,
- Camera SONY HD with remote,
- Video Server,
- Software for Video Server,
- Additional Display 19" on the body of the lamp,
- Additional Display 24" on the body of the lamp,
- Additional Display 42" on the body of the lamp,
- Additional display handle 19",
- Additional display handle 24",
- Additional Display HD 24",
- Wall grip for displays,
- Additional outriggers,
- System for low ceilings,
- Battery for mobile lamp EMALED 300M,
- Charger for mobile lamp EMALED 300M,
- Autoclaved handle for lamp mantle.

