



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Extrication Collar JEROME NecLoc NL200E (red) (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in USA,
- Size: Pediatric (red),
- Ergonomic two-piece clamshell design and retaining strap secure flexion during application,
- Five unique phenotype-driven sizes make it easy to select the correct size,
- Minimal manipulation of patient during application,
- X-ray, MRI & CT transparent,
- Easy to clean,
- Sizes with matching colors for quick recognition,
- Technical Passport.
- Warranty: 30-days

