



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Erlenmeyer Flask VITLAB 500 ml PP no. 9 (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Germany,
- Made from polypropylene,
- Highly transparent,
- With cap,
- Capacity: 1000 ml,
- Graduation: 100 ml,
- Cut NS 45/40,
- Cylindrical neck: 44 mm,
- Blue Scale,
- Optional - sterilisation (autoclaved at 135°).

