



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Beaker Set VITLAB 50 ml PP no. 11 (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Germany,
- Made from polypropylene,
- Highly transparent,
- Convex scale,
- Class B,
- Capacity: 100 ml,
- Graduation: 0.1 ml,
- Height: 200 mm,
- Tolerance: +/- 1 ml,
- Diameter: 24 mm,
- Washable in the temp. up to 20°C.

