



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Pipette POLAMED 0,2 cm³/0,2 ml no. 38 (Ex-Display)

- New,
- Made in Poland,
- Class B,
- Glass,
- Schellbach belt,
- Brown graduation,
- Capacity: 2 cm³ / 2 ml,
- Graduation: 0,002 ml.

