



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Test Tube 10 ml no. 32 (Ex-Display)

- Condition: New: Ex-Display,
- Made in Poland,
- Made from laboratory glass,
- Conical bottom,
- Heat-resistant,
- White Scale,
- Capacity: 10 ml,
- Graduation: 0,1 ml,
- Height: 95 mm,
- Diameter: 12 mm.

