



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Erlenmeyer Flask SIMAX 500 ml no. 15 (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Czech Republic,
- Glass,
- Stable base,
- White scale,
- Capacity: 500 ml,
- Height: 175 mm,
- Neck diameter: 54 mm,
- Flask diameter: 100 mm.

