



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Graduated Flask VITLAB 1000 ml PP (class B) no. 3 (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Germany,
- Made from polypropylene,
- Half-transparent,
- With Cap,
- Capacity: 1000 ml,
- Height: 305 mm,
- Neck GL 24 mm,
- Tolerance: +/- 0,6 ml,
- Wash in the temp. up to 20°C.

