



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Erlenmeyer Flask TERMISIL 200 ml no. 23 (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Poland,
- Made from glass - TERESMIL,
- Narrow neck,
- Capacity: 200 ml,
- Height: 130 mm,
- Neck diameter: 34 mm,
- Flask diameter: 74 mm.

