



Centrifuge DIAMED 1202-23 (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Switzerland,
- Power input: 230 V,
- Frequency: 50 - 60 Hz,
- Power: 70 W,
- Speed: 910 rpm,
- Time: 10 min,
- With break,
- Electro-magnetic lock,
- Automatic balance,
- Very quiet,
- 24 tubes capacity: 5 mm x 70 mm,
- Chamber dimensions: 32 cm,
- Dimensions: 39 x 45 x 26,5 cm,
- Technical Passport,
- Warranty: 30-days

