



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Equipment Trolley 40085 (New)

- New
- Made in Turkey
- Designed for operating and intervention rooms by considering user's requests, made of 304 quality stainless steel
- Two shelves
- Four castors, two of them with brake
- Dimensions: 910 x 560 x 850 mm,
- Technical Passport
- 24 months warranty

