



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Laundry Collection Trolley Uzumcu 40585 / 40590 (New)

- New
- Made in Turkey
- Made of stainless steel
- Detachable and washable fabric collection bags
- Four castors, two of them with brake
- Available in 2 options:
 - Single Bag Model 40585,
 - Double Bag Model 40590,
- Technical Passport
- 24 months warranty

