



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



IV Stand (New)

- New
- Made in Turkey
- Made of stainless steel
- Ergonomic pushing arm for easy transfer
- Accessories basket for storing liquid therapy equipment like serum set, branule
- Infusion devices mountable body structure
- Aluminium base with 5 castors increasing balance
- Hook height can be adjusted between 1200-2100 mm
- Dimensions:
 - Max height: 2100 mm,
 - Base diameter \varnothing 600 mm,
- Technical Passport
- 24 months warranty

