



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Baby Washing Table Uzumcu 40885 (New)

- New
- Made in Turkey
- Made of 304 quality stainless steel designed to be used in new born and infant services by considering user's requests
- One photocell controlled and one manually controlled two separate faucets
- Cabinet with double doors and three drawers
- Baby scale and mattress
- Baby washing sink with dimensions of 500x600x300 mm
- Adjustable 4 legs
- Dimensions: 1650 x 650 x 1100 mm,
- Technical Passport
- 24 months warranty

