



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Wall Mount Kit LMK-02 (New)



With the selective and quality steel, the LMK-02 wall mount assures stability and reliability. It is committed to keep displays firmly on the brick or stud walls, even in harsh environments. The adjustable tilt function provides more comfortable viewing angles. Compatible with various VESA standards, the LMK-02 perfectly suits 32" or larger displays weighing up to 80 kg.

Features:

- Enduring steel provides reliable performance,
- Quick and intuitive installation,
- Flexible tilt adjustment,
- Compatible with VESA standard (Max. 440 x 330 mm),
- Lockable safety bar (padlock not included),
- Weight capacity: 80 kg,
- Technical Passport,
- 24-month warranty

