



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



---

## **RODENSTOCK PRO 1000 Ophthalmic unit (New)**

---

- New,
- Made in Germany,
- Features and advantages of RODENSTOCK PRO 1000 Unit:
  - Robust design,
  - Comfortable, adjustable seat,
  - Timeless, modern design,
  - Column with top LED lighting,
  - Internal RGB LED lighting,
  - Movable table for 2 devices,
  - Mechanical locking of the table,
  - Adjustable phoropter arm,
  - Projector stand,
  - Possible additional cabinet,
- Technical specifications:
  - Rotating and sliding table, for two devices,
  - Locking of the table: mechanical,
  - Max load capacity (table): 45kg (25/20kg),
  - Patient seat: ER 500,
  - Patient seat elevation:  $\pm 20$  cm,
  - Patient seat tilt: up to  $33^\circ$ ,
  - Total weight (without chair): approx. 195/215 kg,
  - Max load capacity (chair): 120 kg,
  - Projector column height: 176 cm/ $\pm 1.5$  cm,
  - The unit is equipped with the following outlets:
    - 2 power sockets (230 V),
    - 1x RS-232, 1x 6-pin (6 V/12 V),
    - 1x LAN,
    - 1x USB 2.0,
  - Reading light: 12 V white LED,
  - Interior light: RGB LED 12 V,
  - Power consumption: 280 VA,
  - Power supply: 110/115 V AC/230 V AC  $\pm 10\%$ ,
- Has a current Technical Passport issued,
- Warranty 12 months,
- Valid technical passport issued,

