



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RODENSTOCK PERISTAT Perimeter (New)

- New,
- Made in Germany,
- Product advantages and features:
 - Automatic measurement,
 - Test validity control,
 - Goldmann III size stimulus,
 - Innovative solutions,
 - Static measurement,
 - Electric chin,
 - Holder for additional lenses,
 - Modern and precise technological solutions,
 - Digital eye tracking,
 - Standardized driver test,
 - Can be configured on almost any PC with Windows operating system,
- Specifications:
 - Measuring bowl (closed bowl with mechanical ventilation):
 - Radius: 30 cm,
 - Measurement range: 50° horizontally (+30° with movable fixation point), 50° vertically,
 - Fixation point: Heijl-Krakau (blind spot method) and CCD camera (digital eye tracking),
 - Chin adjustment: Electrically controlled,
 - User interface: Multilingual,
 - STANDARD TEST FIELDS:
 - Full 50 °,
 - Glare 22 ° / 50 °,
 - Central 22 ° / 30 °,
 - Wide 22 ° / 30 °,
 - Peripheral 30 ° -50 °,
 - Yellow spot 10°,
 - Driving 50 ° / 80 °,
 - SPECIAL STUDY:
 - Bi-Driving: EU license 120 °,
 - Bi-Driving: EU license 160 °,
 - Bi-Driving: full,
 - Industrial medicine,
 - Extended industrial medicine,
 - SOURCE OF BODDER : LED rear projection,
 - SOURCE COLOR : Green,
 - SOURCE SIZE : Goldmann size III (0.43°),
 - BACKLIGHT : 10 asb,
 - DIMENSIONS AND ELECTRICAL REQUIREMENTS:
 - Dimensions (WxDxH): 566 × 633 × 396 mm,
 - Weight: Approximately 18 kg,
 - Power supply: 115 to 230 V AC,
 - Frequency 50/60 Hz,
 - Power consumption: Max 45 VA,
- Has a current Technical Passport issued,
- Warranty 12 months,

