



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical and dental compressor DA7001 (New)

New,
Oil-free,
Power supply 230 V,
Frequency 50 - 60 Hz,
Power 0,75 kW,
Air flow 92 l/min,
Pressure 5 - 7.5 bar,
Tank capacity 30 litres,
Operating sound pressure 51 dB,
Weight 38 kg,
Dimensions: 45 x 45 x 70 cm.
It has a valid Technical Passport,
12-month warranty,
Possibility of financing (Poland only): instalments, leasing, loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

