



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Gynaecological chair DEMETER M01 (New)

New,
Made in Poland,
It has a stable steel construction covered with powder paint,
Additional durability of the chair is ensured by the use of stainless steel in its moving parts,
The bed consists of three parts, including a removable footrest,
DEMETER M01 chair has a fixed height,
The backrest of the armchair is adjustable with a gas spring, after assembling the footrest it is possible to achieve a lying position,
As a standard the armchair is equipped with knee pads adjustable in all planes, a stainless steel container and a folding platform making it easier to get in and out of the armchair,
You can order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available for an extra charge,
Additional options are available (at an extra charge):
Paper hanger,
Caster system,
Hand grips (armrests),
Technical features:
Length with footrest 190, without footrest 139 cm,
Width 65 cm,
Height 85 cm,
Angle of backrest adjustment 00 - 750,
Step height 25 cm,
Weight 70 kg,
Maximum load 150 kg,
It has a current Technical Passport,
Warranty on structure 24 months,
Warranty on hydraulic accessories 12 months,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

