



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical couch AION S02 (New)

New,
Made in Poland,
It has a simple and stable frame made of powder-coated steel profiles,
Fixed height, three-section sleep surface with standard adjustable headrest with cut-out for face and plug,
Headrest is adjustable with a gas spring,
Footrest is adjustable with a rastomat mechanism,
You can order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available at an extra charge,
Additional options are available (at an extra charge):

Hanger for the sleeper,

Castors with brakes,

Bottom shelf,

Armrests,

Technical parameters:

Length 190 cm,

Width 55 or 70 cm,

Height 55 or 70 cm,

empty weight 36 kg,

Maximum load 150 kg,

Headrest angle adjustment -700 to +450,

Adjustment of footrest angle 00 to +350,

It has a valid Technical Passport,

Warranty on structure 24 months,

Warranty on hydraulic accessories 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

