



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical couch ENGO 02 (New)

New,
Made in Poland,
It has a simple and stable frame made of powder-coated steel tubes,
Fixed height, three-section sleep surface with standard adjustable headrest, face hole and plug and adjustable footrest,
Head and footrest are adjustable with a rostomat mechanism,
You can order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available at an extra charge,
Additional options are available (at extra charge):
Hanger for the sleeper,
Castors with brakes,
Bottom shelf,
Technical parameters:
Length 190 cm,
Width 55 or 70 cm,
Height 55 or 70 cm,
Self-supporting weight 34 kg,
Maximum load 150 kg,
Headrest angle adjustment 00 to 450,
Adjustment of footrest angle 00 to 300
It has a valid Technical Passport,
Construction warranty 24 months,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

