



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical couch LETA S02 (New)

New,
Made in Poland,
Has a stable construction of beech wood,
Fixed height, three-section sleep surface with standard adjustable headrest with cut-out for face and plug,
Head and footrest are adjustable using a rostomat mechanism,
It is possible to order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available at an extra charge,
Additional options are available (at extra charge):
Underlay hanger,
Bottom shelf,
Technical parameters:
Length 190 cm,
Width 55 or 70 cm,
Height 55 or 70 cm,
Own weight 32 kg,
Maximum load 150 kg,
Headrest angle adjustment 00 to 450,
Adjustment of footrest angle 00 to 350,
It has a valid Technical Passport,
Warranty: 24 months,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

