



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LETA R02 medical couch (New)

New,
Made in Poland,
It has a stable construction of beech wood,
Adjustable height and three-section sleep surface with standard adjustable headrest with cut-out for face and plug,
Height adjustment of the couch is stepped,
The headrest and the middle and back parts of the couch are adjustable using a rastomat mechanism,
The upholstery can be ordered free of charge in any colour you choose, except for "PRESTIGE" colours which are available for an extra charge,
Additional options are available (at extra charge):
Underlay hanger,
Bottom shelf,
Technical parameters:
Length 190 cm,
Width 55 or 70 cm,
Height 55 to 70 cm,
Weight 35 kg,
Maximum load 150 kg,
Headrest angle adjustment 00 to 450,
Adjustment of footrest angle 00 to 350,
It has a valid Technical Passport,
Warranty: 24 months,
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

