



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BLOODY 01 blood collection chair (New)

New,

Made in Poland,

The design makes it easy to carry out the injection procedure on the left or right forearm of the person subject to the procedure,

It is used in doctor's surgeries, health clinics, hospitals and laboratories,

The chair has a stable, powder-coated steel construction, adjustable backrest and two armrests adjustable in all planes,

The armchair can be fitted with a headrest for an additional charge,

It is possible to order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available at an extra charge,

Technical parameters:

Height 100 cm,

Width 50 cm,

Depth 80 cm,

Seat height 50 cm,

Maximum load 150 kg,

own weight 24 kg,

Headrest height adjustment 15 cm

Height adjustment of armrests 25 cm,

It has a current Technical Passport,

Warranty: 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

