



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Stool ERIS 04 (New)

New,

Polish production,

It has a stable aluminium base, comfortable upholstered seat and practical footrest,

Height adjustment is done with a gas spring,

Height adjustment is done by gas spring. Backrest available in two variants, wide or narrow (20x10,20x20)

Stool is available on wheels or feet,

It is possible to order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available at an extra charge,

Technical parameters:

Height on feet 59 - 84 cm,

Height on wheels 64 - 89 cm,

Base ø 65 cm,

Seat ø 35 cm,

12-month warranty,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

