



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Endoscopic Camera System AESCULAP PV 206 (Reconditioned)

- Condition: Reconditioned - Very Good,
- Made in Germany,
- Power input: 230 V,
- Frequency: 50 - 60 Hz,
- Power: 35 VA,
- Adjustable with rigid endoscopes,
- Sealed control panel,
- Video Signal Outputs:
 - S-VIDEO,
 - VIDEO,
- Dimensions: 30,5 x 30,5 x 8,5 cm,
- Technical Passport,
- 30-days warranty

