



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Baby examination table GAJA 1 (New)

New,

Made in Poland,

It is designed to perform examinations and care of infants in clinics and doctor's offices,

The infant examination table is made of powder-coated steel pipes.

It has a fixed height,

The top is protected on three sides with high foam and at the bottom there is a comfortable removable mattress, which makes it easier to keep the table clean.

Additional options are available (at an extra charge):

Bottom shelf

Side shelf

Drawer

Cabinet

Hanger for the underlay,

Technical parameters:

Total width of the top 80 or 100 cm,

Total depth of the top 65 or 75,

Total height 79 or 89 cm,

Weight 30 kg,

Maximum load 50 kg,

Maximum load of the side shelf 25 kg,

It has a valid technical passport

Construction warranty 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

