



Transfer Chair NITROCARE NTS X6 (New)

- New
- Made in Turkey
- Easy of use with ergonomic design, light structure
- Excellent maneuverability
- Adjustable arm supports
- Armrests with gradual opening and closing
- Central brake bar
- Adjustable headrest, backrest
- Polyurethane backrest and seating
- Dimensions;
 - Height: 1160 mm
 - Depth: 1000 mm,
 - Width: 643 mm,
 - Seat depth: 440 mm,
 - Backrest width: 500 mm,
 - Armrests: 355 x 425 x 65 mm,
- Technical Passport
- 24 months warranty

