



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Laryngology chair PROMAT NG (New)

New,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Standard colour is WHITE/Black,
Functionality of the chair is ensured by removable multifunctional headrest,
Possibility of installing a holder for mounting a paper roll,
Equipped with 4 electric motors,
Controlled by remote control or foot control:
Height adjustment in the range: 615 - 890 mm,
Backrest adjustment from: 90° - 180°,
Footrest adjustment from: 180° - 255°,
Seat adjustment from: -2° - 32°,
Manual adjustment:
Headrest adjustment: 180°,
Adjustment of footrest extension: 180°,
Adjustment of armrests: 180°,
Maximum load capacity of the chair: 240 kg,
Retractable castors for easy manoeuvrability,
Additional equipment at extra charge:
Holder for mounting a paper roll,
Foot control,
Accessory holder,
IV holder with four hooks,
Seat rail system,
Back rail system,
Ability to change the colour of the upholstery,
Noise level: 60 dBA,
Weight: 90 kg,
Dimensions:
Headrest: 274 x 415 mm,
Backrest: 710 x 600 mm,
Seat: 510 x 600 mm,
Footrest: 305 x 530 mm,
Armrest: 555 x 95 mm,
Footrest extension: 170 x 450 mm,
Overall length from 1950 to 2230 mm,
Overall width: 600/845 mm (including armrests),
With current Technical Passport,
Warranty: 24 months,
Finance options (Poland only): Instalments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

