



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OTODYNAMICS ECHOCHECK (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in UK,
- Power input: 230 V or battery 2,5 AH,
- Frequency range: 1,6 kHz do 2,8 kHz,
- Signal range: -6 dB SPL to 55 dB SPL,
- Automatic tests result,
- LED Display,
- Compatible with PC,
- Comes with:
 - Device,
 - Probe,
 - Charger,
 - CD software,
 - Case,
- Weight: 370 g,
- Dimensions: 19 x 9 x 3 cm,
- Technical Passport,
- Warranty: 30-days

