



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Ophthalmic unit Classic Comfort PLUS (New)**

New,

Polish production,

Characteristics:

The unit consists of table top for 2 devices, a base, a place for the projector and an electrically raised seat

It has a high quality actuator made by Swiss company Schumo,

The seat is upholstered with skaita in any colour from the pattern,

It can be locked in any position,

The unit is equipped with lighting in the form of a halogen lamp with an option to regulate light intensity and optionally: a balance phoropter arm,

Standard equipment:

Footrest,

Lamp,

Electrically raised seat,

Electronic control panel,

Optional equipment (at extra charge):

Phoropter arm,

Technical data:

Supply voltage 230 V, 50 Hz

Power consumption 460 W

Voltage available on tabletop 230 V, 50 Hz + option 6-12 V

Maximum load 80 kg

Weight 180 kg

Height 184 cm

Dimensions 1240 x 1594 cm

It has a current Technical Passport,

12 months warranty,

Finance options (Poland only): instalments, leasing, loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

