



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Ophthalmic unit Classic Comfort (New)

New,

Polish production,

Characteristics:

Ergonomic Schumo power-lifting armchair, equipped with footrest, armrests upholstery made of leather in a colour of your choice,

Elegant base,

Stable base,

Intuitive operation with a handy control panel,

Footrest provides comfortable support for the patient's feet,

Modern design and high quality materials ensure long term satisfactory use and easy cleaning,

The presented unit consists of a sliding table top for 2 devices, a base, a place for a projector and an electrically raised chair,

Standard equipment:

Large top electromagnetically locked in the desired position,

Power sockets: 230V and 6V / 12V for slit lamp,

Shelf under the projector,

Halogen lamp with light intensity regulation

Additional equipment (at extra charge):

Possibility of installing a balance phoropter arm and an intuitive control panel for the COMF ophthalmic unit,

Armrests,

Technical data:

Supply voltage 230 V, 50 Hz,

Power consumption 460 W,

Voltage available on tabletop 230 V, 50 Hz + option 6-12 V,

Maximum load 80 kg,

Weight 180 kg,

Height 184 cm,

Dimensions 1240 x 1594 cm,

Warranty 12 months,

Financing option (Poland only): instalments, leasing, loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE



