



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Ophthalmologic Table SE-1 (New)

- New,
- Made in Poland,
- Compatible with ophthalmologic diagnostic devices,
- Stable base with aluminium feets - electric actuator from SCHUMO (Switzerland) with high range of height adjustment,
- Countertop 600 x 400 mm (length for 1 devices),
- Socket 230V mounted under the countertop,
- Countertop color: Standard Grey / Optional: colors from RAL palette,
- Optional: on casters + additional shelf for the keyboard,
- Technical specification:
 - Countertop: 600 x 440 mm,
 - Height: Min. 700, max. 900 mm,
 - Weight: 30 kg,
 - Power input: 230 V, 50 - 60 Hz,
 - Maximum load: 60 kg,
 - Classification: Medical device I Class type B,
- Technical Passport,
- 12-month warranty

