



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LID WOOD "Blacklight" Fluorescent light (New)

- New,
- Made in France,
- LID - WOOD "Blacklight" lamp is a diagnostic tool that emits ultraviolet light used in many fields,
- In dermatology, it is used to examine the skin of patients to detect several types of problems,
- It is particularly useful in the diagnosis of erythema, skin fungal infections, pigmented lesions, bacterial skin infections, as well as pigmentation problems (such as acquired vitiligo), melanomas, dermatoses, etc.,
- This portable lamp is lightweight and easy to use. The power socket is equipped with an electronic ballast, and the bulb is protected by a synthetic glass shade.
- Technical specifications:
 - Wavelength: 370 nm,
 - Fluorescent lamp 9 W (300200PH),
 - Weight: 1 kg,
 - Lamp power - 9 W
 - Life span of the light tube minimum 3,000 hours of operation,
 - Power cord length: 280 cm,
 - Available colors: black,
 - Dimensions: 34 cm x 10 cm
- Has a current Technical Passport issued,
- Warranty: 3 years,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

